

# Protein expression for crystal generation

 Sean William Fanning

Updated date: Oct 12, 2020

 An abbreviated version of this protocol was published in eLIFE in Feb 2016

Estrogen receptor alpha somatic mutations Y537S and D538G confer breast cancer endocrine resistance by stabilizing the activating function-2 binding conformation

DOI: 10.7554/eLife.12792

## Detailed protocol

## Updated Expression and Purification Protocol

Running/Lysis Buffer 1L keep refrigerated.

- 50 mM HEPES pH 8 (Tris pH 8 is fine too)
- 500 mM NaCl
- 0.5 mM TCEP
- 20 mM Imidazole
- Fill with DI H<sub>2</sub>O to 900 mL
- Add 100 mL glycerol while stirring.

Elution Buffer 500 mL

- 50 mM HEPES pH 8
- 1 M NaCl
- 500 mM TCEP
- 500 mM Imidazole
- Fill with DI H<sub>2</sub>O to 450 mL
- Add 50 mL glycerol (under sink) while stirring.

Resuspension

- 10 mL cold running buffer per gram cell paste.
- 1 EDTA-Free pic tab per 50 mL lysis (Pierce Cat#A32965).
- Thoroughly resuspend on ice by pipetting.

Lysis with Sonication

- Per 50 mL lysate, 5 s on/off cycles for 5 minutes of total sonicating at 80% power on ice. Sonicator used is a Fisherbrand Model 120 Sonic Dismembrator with a horn suitable for 50 mL sonication volumes.

Clarifying Lysate

- Centrifuge at >18,000 xg for 30 minutes at 4°C using the fixed angle rotor. Lysate should be clear and yellow.
- Transfer to a new tube to separate soluble and insoluble fractions, make sure that no insoluble debris transfers over.

Purification

- Purification was performed using a BioRad NGC Quest FPLC with a sample pump.
- We use a BioRad Nuvia 5 mL IMAC column for affinity purification. FPLC should have the running/lysis buffer loaded on pump A and the elution buffer on pump B.
- Clean the column by flushing 3 CV of 100% pump B.
- Equilibrate the column with 5 CV of pump A or until UV baseline is reached.
- Place sample pump line into the clarified lysate and pull 5 mL through the pump to prime it and purge any air.
- Load the sample directly onto the column, wash with 100 pump A until baseline is reached, do a linear gradient elution from 0 to 100% pump B.
- Pick peak fractions and run 20  $\mu$ L (with 20  $\mu$ L load dye) on a SDS-PAGE gel along with ladder.
- Combine fractions with the correct band corresponding to protein of interest take a reading on the nanodrop. A<sub>260</sub>/A<sub>280</sub> should be <0.8. If greater, add DNase. Use extinction coefficient to determine molarity.

TEV Cleavage and Dialysis

- Add 1:5,000 HisTev, place in dialysis membrane, and dialyze against 4 L of running buffer overnight at 4°C with stirring.

Reverse Ni-NTA

- Load dialyzed protein onto the IMAC column as before, collect flow through. Run 20 mL running buffer over column and collect. In a separate fraction, run 100% pump B/elution buffer to elute the cleaved His tag and His-TEV.
- While this is running make SEC buffer, which is running buffer without imidazole. Degas for at least 30 minutes with stirring. Once reverse Ni-NTA protocol is complete, start equilibrating the SEC column at 1 mL/min for at least 1 CV.

Gel Filtration

- Using a spin concentrator with a 30 kDa molecular weight cutoff, concentrate the protein to at least 5 mg/mL using the filtration concentrators. Use the swinging bucket rotor at 4°C, 4,000xg for 5-15 minutes at a time. Concentrated to no less than 1 mL, ideally to 5 mL. Protein should not be concentrated to more than 20 mg/mL, make sure you are taking nanodrop readings throughout the concentrating. Ideally, the A<sub>260</sub>/A<sub>280</sub> ratio will decrease as the A<sub>280</sub> increases. Remember that the absorbance isn't always 1:1 with the A<sub>280</sub> given the extinction coefficient of the protein, but in general you need an A<sub>280</sub> of at least 3 to see fractions on the SEC.
- Run fractions on the gel. ER LBD typically elutes around 30-60 kDa on the SEC. Concentrate to 5-20 mg/mL and flash freeze for later use. If time, set up crystal trays right away.

**How to cite:** (Readers should cite both the Bio-protocol preprint and the original research article where this protocol was used)

1. Fanning, S. W. (2020). Protein expression for crystal generation. Bio-protocol Preprint. [bio-protocol.org/prep538](https://bio-protocol.org/prep538).
2. Fanning, S. W., Mayne, C. G., Dharmarajan, V., Carlson, K. E., Martin, T. A., Novick, S. J., Toy, W., Green, B., Panchamukhi, S., Katzenellenbogen, B. S., Tajkhorshid, E., Griffin, P. R., Shen, Y., Chandarlapaty, S., Katzenellenbogen, J. A. and Greene, G. L. (2016). Estrogen receptor alpha somatic mutations Y537S and D538G confer breast cancer endocrine resistance by stabilizing the activating function-2 binding conformation. eLIFE. DOI: [10.7554/eLife.12792](https://doi.org/10.7554/eLife.12792)

